

AF/LE Purified Anti-Mouse CD45R/B220 Antibody[RA3.3A 1/6.1]

Catalog No. E-AB-F11120

Note: Centrifuge before opening to ensure complete recovery of vial contents.**Description**

Reactivity	Mouse
Host	Rat
Isotype	Rat IgM, κ
Clone	RA3.3A 1/6.1
Conjugation	Unconjugated
Buffer	0.2 μ m filtered in PBS, pH 7.2. Azide Free (AF)/Low Endotoxin (LE): Contains no stabilizers or stabilizers. Endotoxin level is < 2 EU/mg as Determined by LAL gel clotting assay.

Preparation & Storage**Storage** This product can be stored at 2-8°C for 12 months. Do not freeze.**Background**

CD45R, also known as B220, is an isoform of CD45. It is a member of the protein tyrosine phosphatase (PTP) family with a molecular weight of approximately 180-240 kD. CD45R is expressed on B cells (at all developmental stages from pro-B cells through mature B cells), activated B cells, and subsets of T and NK cells. CD45R (B220) is also expressed on a subset of abnormal T cells involved in the pathogenesis of systemic autoimmunity in MRL-Fas^{lpr} and MRL-Fas^{gld} mice. It plays a critical role in TCR and BCR signaling. The primary ligands for CD45 are galectin-1, CD2, CD3, and CD4. CD45R is commonly used as a pan-B cell marker; however, CD19 may be more appropriate for B cell specificity.

For Research Use Only